

NMS

- Light weight (Transformerless)
- With heat filter

Digital Cold Light Source

(24V x 250 watts twin bulbs)

Digital Cold Light Source



FEATURES

- Light weight, Metal casing.
- Light intensity is digitally controlled by a microprocessor based device.
- Inbuilt heat filter to provide genuine cold light.
- Due to Genuine cold light, casing and light adapter remain cold during operation.
- Easy lamp change.
- Adheres to high international safety standard.

TECHNICAL SPECIFICATION (With SMPS)

Type of Lamp	Halogen reflector lamp 24V - 250 W
Electrical Supply	160V to 270V AC, 50Hz
Power consumption	350VA
Fuse	3A
Dimensions	32cm(d) x 25cm(w) x 13cm(h)
Weight	5 kg.
Product Code	NMS-LH-SM

TECHNICAL SPECIFICATION (With Transformer)

Type of Lamp	Halogen reflector lamp 24V - 250 W
Electrical Supply	230V AC \pm 10% 50Hz
Power consumption	350VA
Fuse	5A
Dimensions	32cm(d) x 25cm(w) x 13cm(h)
Weight	9 kg.
Product Code	NMS-LH-TR