

NIDHI⁺

For mini PCNL & Ureteroscopy

Easyflow Pump

Easyflow Pump



FEATURES

- This classical Easyflow pump is a unique device for urological applications i.e.
- The flow and pressure can be selected from front panel.
- By the use of provided silicone tubing set precise pre selected pressure and flow can be transmitted to surgical device.
- Over pressure warning alarm and cut-off facility to provide high level safety for patient and surgeon.
- Tubing set silicone - autoclavable.

Application : It is used to dilate ureter during ureteroscopy. It is also used for mini PCNL procedure .

TECHNICAL SPECIFICATION

Power supply	160 to 270V AC, 50Hz
Fuse	2 Amps
Flow Rate	0 - 1200 ml /min.
Pressure Range	0 - 400 mm Hg.
Temperature Range	10° to 40°
Dimension	26cm(d) x 36cm(w) x 17cm(h)
Weight	6.5 kg. max.