


UROLOGY DISPOSABLES



UROLOGY DISPOSABLES

Nidhi Meditech System Provides genuine endo urological disposables like Double Pigtail Ureteral stent, Mono 'J' Stent, Endopylotomy stent, Ureteric catheters etc. All the endourological disposables manufactured as per ISO 13485 : 2003 and CE 1023 standards.

All the products manufactured from tecoflex grade Polyurethanes according to ISO-10993 series


In general, the ease of processing and high strength of polyurethane makes them the material of choice in medical devices. Polyurethanes pass the tests of biocompatibility and chemical stability. They can be easily loaded with radiopaque material for X-ray detection. However, the major advantages of polyurethane are their thermolability that makes them very soft inside the body, minimizing the sense of discomfort associated with implantable foreign body. Nidhi double loop stent is made of tecoflex, which has the maximum degree of softening inside the body. It has a wide margin of durometer extending from very soft to very hard with different sizes and length according to the physician specification.

DOUBLE PIGTAIL URETERAL STENT

FEATURES


1. Made of medical grade biocompatible techoflex polyurethane for Smooth Passage radio-opaque implant.
2. Durable with superior bio-compatibility
3. Thin wall (provides potency and good drainage).
4. Atraumatic insertion is revealed quite smooth.
5. Minimized tissue irritation with easy but rapid atraumatic insertion.
6. Cost effective.

DOUBLE PIGTAIL URETERAL STENT WITH PUSHER



Size	Product No.	Length	Type	Compatible Guide Wire..
4.5 FR	NMS-DJ-4.5 FR	26 cms.	(OEW/CEW)	0.025"
4.8 FR	NMS-DJ-4.8 FR	26 cms.	(OEW/CEW)	0.032"
5 FR	NMS-DJ- 5 FR	26 cms.	(OEW/CEW)	0.035"
5.5 FR	NMS-DJ- 5.5 FR	26 cms.	(OEW/CEW)	0.035"
6 FR	NMS-DJ- 6 FR	26 cms.	(OEW/CEW)	0.038"


DOUBLE PIGTAIL URETERAL STENT WITH PUSHER & SUTURE



Size	Product No.	Length	Type	Compatible Guide Wire..
4.5 FR	NMS-DJ-ST-4.5 FR	26 cms.	(OEW/CEW)	0.025"
4.8 FR	NMS-DJ-ST-4.8 FR	26 cms.	(OEW/CEW)	0.032"
5 FR	NMS-DJ-ST- 5 FR	26 cms.	(OEW/CEW)	0.035"
5.5 FR	NMS-DJ-ST- 5.5 FR	26 cms.	(OEW/CEW)	0.035"
6 FR	NMS-DJ-ST- 6 FR	26 cms.	(OEW/CEW)	0.038"

UROLOGY DISPOSABLES


BIO FLEX DOUBLE PIGTAIL URETERAL STENT WITH PUSHER



Application	Applied for internal drainage from the ureteropelvic junction to the bladder, as a provisional measure.
Features	Made of medical grade biocompatible soft techoflex polyurethane for Smooth Passage radio-opaque implant.
Use	For single use
Packing	Tyvek pouch as per EN-868-5-1999 standard
Sterilizations	With Ethylene Oxide EN-550 series standard

Size	Product No.	Length	Type	Compatible Guide Wire..
5 FR	NMS-BFDJ- 5 FR	26 cms.	(OEW/CEW)	0.035"
6 FR	NMS-BFDJ- 6 FR	26 cms.	(OEW/CEW)	0.038"

BLACK SILICON DOUBLE PIGTAIL URETERAL STENT WITH PUSHER



Size	Product No.	Length	Type	Compatible Guide Wire..
5 FR	NMS-BSDJ- 5 FR	26 cms.	(OEW/CEW)	0.035"
6 FR	NMS-BSDJ- 6 FR	26 cms.	(OEW/CEW)	0.038"

Application	Applied for internal drainage from the ureteropelvic junction to the bladder, as a provisional measure.
Features	Made of medical grade biocompatible silicon material for Smooth Passage radio-opaque implant.
Use	For single use
Packing	Tyvek pouch as per EN-868-5-1999 standard
Sterilizations	With Ethylene Oxide EN-550 series standard

- * Please specify in your order product code, Single Loop / Multi Loop
- * Any special size & length can be available on request.
- * OEW = Both End Open White, CEW = One End Close White.

Product Code: Product No. + Length + Type

UROLOGY DISPOSABLES

STONE RETRIEVAL BASKET - HELICAL (S.S.)

STONE RETRIEVAL BASKET


- Application** To remove multiple stone fragments from ureter with semi rigid endoscope.
- Features** Maximum patient safety and compatibility. easy to operate.
- Use** For single use
- Packing** Blister tray with tyvek paper as per EN-868-5-1999 standard
- Sterilizations** With Ethylene Oxide EN-550 series standard



Size (In FR)	Product Code	Length
3 / 4 WIRE	NMS-SS-H-3FR-4FR-70/90	70/90 cms.
4 / 4 WIRE	NMS-SS-H-4FR-4FR-70/90	70/90 cms.
5 / 4 WIRE	NMS-SS-H-5FR-4FR-70/90	70/90 cms.

STONE RETRIEVAL BASKET - SEGURA (S.S.)

- Application** To remove multiple stone fragments from ureter with semi rigid endoscope.
- Features** Maximum patient safety and compatibility. Easy to operate.
- Use** For single use
- Packing** Blister tray with tyvek paper as per EN-868-5-1999 standard
- Sterilizations** With Ethylene Oxide EN-550 series standard



Size (In FR)	Product Code	Length
3 / 4 WIRE	NMS-SS-H-3FR-4FR-"0"-70/90	70/90 cms.
4 / 4 WIRE	NMS-SS-H-4FR-4FR-"0"-70/90	70/90 cms.
5 / 4 WIRE	NMS-SS-H-5FR-4FR-"0"-70/90	70/90 cms.

Size (In FR)	Product Code	Length
3 / 4 WIRE	NMS-SS-S-3FR-4FR-"0"-70/90	70/90 cms.
4 / 4 WIRE	NMS-SS-S-4FR-4FR-"0"-70/90	70/90 cms.
5 / 4 WIRE	NMS-SS-S-5FR-4FR-"0"-70/90	70/90 cms.

H=Helical

S=Segura


"0" = '0' TIP

UROLOGY DISPOSABLES

STONE RETRIEVAL BASKET - HELICAL (NITINOL)



Application To remove multiple stone fragments from ureter with semi rigid/flexible endoscope.
Features Maximum patient safety and compatibility. easy to operate.
Use For single use
Packing Blister tray with tyvek paper as per EN-868-5-1999 standard
Sterilizations With Ethylene Oxide EN-550 series standard




Size (In FR)	Product Code	Length
3 / 4 WIRE	NMS-NT-H-3FR-4FR-70/90	70/90 cms.
4 / 4 WIRE	NMS-NT-H-4FR-4FR-70/90	70/90 cms.
5 / 4 WIRE	NMS-NT-H-5FR-4FR-70/90	70/90 cms.

STONE RETRIEVAL BASKET - SEGURA (NITINOL)



Application It is used for removing the large stones from ureter with semi rigid/flexible endoscope.
Features Retrieving the large stones, smoothly.
Use For single use
Packing Blister tray with tyvek paper as per EN-868-5-1999 standard
Sterilizations With Ethylene Oxide EN-550 series standard



Size (In FR)	Product Code	Length
3 / 4 WIRE	NMS-NT-H-3FR-4FR-"0"-70/90	70/90 cms.
4 / 4 WIRE	NMS-NT-H-4FR-4FR-"0"-70/90	70/90 cms.
5 / 4 WIRE	NMS-NT-H-5FR-4FR-"0"-70/90	70/90 cms.

Size (In FR)	Product Code	Length
3 / 4 WIRE	NMS-NT-S-3FR-4FR-"0"-70/90	70/90 cms.
4 / 4 WIRE	NMS-NT-S-4FR-4FR-"0"-70/90	70/90 cms.
5 / 4 WIRE	NMS-NT-S-5FR-4FR-"0"-70/90	70/90 cms.

UROLOGY DISPOSABLES

MALECOT NEPHROSTOMY CATHETER



Size	Length	Product Code
10 fr	30 cms	NMS-MC-10FR-30
12 fr	30 cms	NMS-MC-12FR-30
14 fr	30 cms	NMS-MC-14FR-30
16 fr	30 cms	NMS-MC-16FR-30
18 fr	30 cms	NMS-MC-18FR-30


*** Any special size & length can be available on request.**

MALECOT NEPHROSTOMY CATHETER (WITH NEEDLE)

- Application** The Malecot nephrostomy catheter is placed percutaneously through an existing nephrostomy tract. It is used for nephrostomy drainage in the renal pelvis.
- Features** Smooth passage in tortuous ureter.
- Use** For single use
- Packing** Tyvek pouch as per EN-868-5-1999 standard
- Sterilizations** With Ethylene Oxide EN-550 series standard


Size	Length	Product Code
10 fr	30 cms	NMS-MC-N-10FR-30
12 fr	30 cms	NMS-MC-N-12FR-30
14 fr	30 cms	NMS-MC-N-14FR-30
16 fr	30 cms	NMS-MC-N-16FR-30
18 fr	30 cms	NMS-MC-N-18FR-30

*** Any special size & length can be available on request.**



ENDOPYLOTOMY STENT

Product Code : ENDOPYLOTOMY STENT




- Application** Used for Stenting the ureteropelvic junction & ureter following a percutaneous endopyelotomy while providing both nephrostomy & ureteral drainage
- Use** For single use
- Packing** Tyvek pouch as per EN-868-5-1999 standard
- Sterilizations** With Ethylene Oxide EN-550 series standard

5 fr / 10 fr	26 cms	NMS-EPS-5/10FR26
6 fr / 12 fr	26 cms	NMS-EPS-6/12FR26
7 fr / 14 fr	26 cms	NMS-EPS-7/14FR26

UROLOGY DISPOSABLES

DOUBLE PIGTAIL URETERAL STENT SET (WITH PUSHER, GUIDE WIRE & CLAMP)



Application	Applied for internal drainage from the ureteropelvic junction to the bladder, as a provisional measure.
Features	Made of soft telecoflex polyurethane for Smooth Passage radio-opaque implant.
Contents	Pusher, Suture, Clamp with PTFE coated Guidewire
Use	For single use
Packing	Tyvek pouch as per EN-868-5-1999 standard
Sterilizations	With Ethylene Oxide EN-550 series standard

Size	Product No.	Length	Type	Compatible Guide Wire..
4.5 FR	NMS-DJ-SET-4.5 FR	26 cms.	(OEW/CEW)	0.025"
4.8 FR	NMS-DJ-SET-4.8 FR	26 cms.	(OEW/CEW)	0.032"
5 FR	NMS-DJ-SET- 5 FR	26 cms.	(OEW/CEW)	0.035"
5.5 FR	NMS-DJ-SET- 5.5 FR	26 cms.	(OEW/CEW)	0.035"
6 FR	NMS-DJ-SET- 6 FR	26 cms.	(OEW/CEW)	0.038"

Product Code: **Product No.** + **Length** + **Type**

- * Please specify in your order product code, Single Loop / Multi Loop
- * Any special size & length can be available on request.
- * OEW = Both End Open White, CEW = One End Close White.

URETERIC CATHETER

URETERIC CATHETER (OPEN END)

Application	Used for drainage as well as retrograde pyelogram by open end tip it is placed over a guide wire for negotiation to tortuous ureter.
Features	Smooth passage in tortuous ureter.
Use	For single use
Packing	Tyvek pouch as per EN-868-5-1999 standard
Sterilizations	With Ethylene Oxide EN-550 series standard

Size	Product Code	Length
3 FR	NMS-UC-3 FR 70 (OEG)	70 cms.
4 FR	NMS-UC-4 FR 70(OEG)	70 cms.
5 FR	NMS-UC-5 FR 70 (OEG)	70 cms.
6 FR	NMS-UC-6 FR 70(OEG)	70 cms.

URETERIC CATHETER (CLOSED END)

Application	Used for drainage and retrograde pyelogram.
Features	Best in its inherent quality providing drainage and retrograde pyelogram.
Use	For single use
Packing	Tyvek pouch as per EN-868-5-1999 standard.
Sterilizations	With Ethylene Oxide EN-550 series standard

Size	Product Code	Length
3 FR	NMS-UC-3 FR 70 (CEG)	70 cms.
4 FR	NMS-UC-4 FR 70 (CEG)	70 cms.
5 FR	NMS-UC-5 FR 70 (CEG)	70 cms.
6 FR	NMS-UC-6 FR 70 (CEG)	70 cms.



CEG = One End Close Green, **OEG** = Both End Open Green

* Any special size & length can be available on request.

UROLOGY DISPOSABLES

Product Code : NMS-PCN

MONO-J-STENT



- Application** By direct puncture, it provides temporary drainage of urine, from the kidney. This is done in presence of distal obstruction.
- Features** Provides drainage from kidney, smoothly.
- Use** For single use
- Packing** Tyvek pouch as per EN-868-5-1999 standard
- Sterilizations** With Ethylene Oxide EN-550 series standard

MONO-J-STENT


Size	Product Code	Length
6 FR	NMS-PCN-6 FR- 22	22 cms.
6 FR	NMS-PCN-6 FR 30	30 cms.
8 FR	NMS-PCN-8 FR 22	22 cms.
8 FR	NMS-PCN-8 FR - 30	30 cms.
9 FR	NMS-PCN-9 FR 22	22 cms.
9 FR	NMS-PCN-9 FR 30	30 cms.
10 FR	NMS-PCN-10FR 22	22 cms.
10 FR	NMS-PCN-10FR 30	30 cms.
12 FR	NMS-PCN-12FR 22	22 cms.
12 FR	NMS-PCN-12FR 30	30 cms.
14 FR	NMS-PCN-14FR 22	22 cms.
14 FR	NMS-PCN-14FR 30	30 cms.

MONO-J-STENT WITH TROCAR



Size	Product Code	Length
6 FR	NMS-PCN-N-6 FR- 22	22 cms.
6 FR	NMS-PCN-N-6 FR 30	30 cms.
8 FR	NMS-PCN-N-8 FR 22	22 cms.
8 FR	NMS-PCN-N-8 FR - 30	30 cms.
9 FR	NMS-PCN-N-9 FR 22	22 cms.
9 FR	NMS-PCN-N-9 FR 30	30 cms.
10 FR	NMS-PCN-N-10FR 22	22 cms.
10 FR	NMS-PCN-N-10FR 30	22 cms.
12 FR	NMS-PCN-N-12FR 22	22 cms.
12 FR	NMS-PCN-N-12FR 30	30 cms.
14 FR	NMS-PCN-N-14FR 22	22 cms.
14 FR	NMS-PCN-N-14FR 30	30 cms.

MONO-J-STENT KIT (WITH NEEDLE, GUIDE WIRE, DILATOR SET & CLAMP)

Product Code : NMS-PCN-KIT



- Application** By direct puncture, it provides temporary drainage of urine, from the kidney. This is done in presence of distal obstruction.
- Features** Provides drainage from kidney, smoothly.
- Contents** Mono-J-Stent, Trocar, 6 Dilators, Guide Wire, Clamp, Scalpel
- Use** For single use
- Packing** Tyvek pouch as per EN-868-5-1999 standard
- Sterilizations** With Ethylene Oxide EN-550 series standard



UROLOGY DISPOSABLES

Application By direct puncture, it provides temporary drainage of urine, from the kidney. This is done in presence of distal obstruction.

Features Provides drainage from kidney, smoothly.

Contents Mono-J-Stent, Trocar, Guide Wire, Clamp, Scalpel

Use For single use

Packing Tyvek pouch as per EN-868-5-1999 standard

Sterilizations With Ethylene Oxide EN-550 series standard


MONO-J-STENT SET (WITH NEEDLE, GUIDE WIRE & CLAMP)

Product Code : NMS-PCN-SET



Product Code : NMS-FDS

FASCIAL DILATORS SET (PCN)



Application Used for fascial dilation prior to percutaneous nephrolithotomy and for stone manipulation.

Features Smooth surfaces for ease of introduction. Uniform taper to reduce trauma.

Size 6 FR - 16 FR

Length 22 cm

Use For single use

Packing Tyvek pouch as per EN-868-5- series standard


Sterilizations With Ethylene Oxide EN-550 series standard

UROLOGY DISPOSABLES

Product Code : NMS-UDS-70

URETERAL DIALATORS SET


Application	Used to dilate the intramural tunnel and ureter in one step. May also be used to inject contrast medium to determine sub mucosal tunneling or anesthetic jelly to aid in dislodging ureteral stones.
Features	-Smooth surfaces for ease of introduction. -Uniform taper to reduce trauma. -Use with 0.038" guide wire
Size	6 FR - 16 FR
Length	70 cm
Use	For single use
Packing	Tyvek pouch as per EN-868-5-series standard
Sterilizations	With Ethylene Oxide EN-550 series standard



Product Code : NMS-UDS-40

URETHRAL DIALATORS SET

Application	Used for dilation of urethra prior to ureterorenoscopy and also for sequential (step by step) dilation of the intramural tunnel and urethra and / or stone manipulation.
Features	Smooth surfaces for ease of introduction. Uniform taper to reduce trauma.
Size	6 FR - 16 FR / 6 FR - 24 FR
Length	40 cm
Use	For single use
Packing	Blister tray with tyvek paper as per EN-868-5-series standard
Sterilizations	With Ethylene Oxide EN-550 series standard



UROLOGY DISPOSABLES

NMS-ST

STONE TRAP



- Application It helps increase procedural efficiency by preventing retro-pulsion of stones up the ureter and into the kidney .
- Features to prevent stone migration in the kidney during procedure.
- Use For single use
- Packing Tyvek pouch as per EN-868-5- series standard
- Sterilizations With Ethylene Oxide EN-550 series standard

NMS-UAS

URETHRAL ACCESS SHEATH



- Application used to facilitate flexible ureteroscopy in the treatment of urolithiasis.
- Features reducing operative times, allowing multiple passes of instruments, and protecting the ureteroscope during flexible ureteroscopy.
- Use For single use
- Size 9.5 / 11.5 FR
- Length 35 cm / 45 cm
- Packing Tyvek pouch as per EN-868-5- series standard
- Sterilizations With Ethylene Oxide EN-550 series standard

NMS-ARDS

AMPLATZ RENAL DILATOR SET



- Application Used for stretching or enlarging the fascial tissue prior to a nephrostomy procedure and to introduce the Amplatz sheath which is intended for percutaneous dilation of tissue surrounding a dilator which is previously placed
- Features radiopaque and has a smooth surface.
- Use For single use
- Packing Tyvek pouch as per EN-868-5- series standard
- Sterilizations With Ethylene Oxide EN-550 series standard

UROLOGY DISPOSABLES

NMS-AMP



AMPLATZ SHEATH

Application	used for smooth passage of surgical instruments into the nephrostomy tract.
Features	smooth surface.
Size	18fr/20fr/22fr/24fr/26fr/28fr/30/
Length	17 cms
Use	For single use
Packing	Tyvek pouch as per EN-868-5- series standard
Sterilizations	With Ethylene Oxide EN-550 series standard

NMS-IPN



I.P. NEEDLE

Application	Intended to make the initial percutaneous puncture to access the kidney and to enable the introduction of a guidewire during urological procedures (PCN, PCNL, Mini PERC & Micro PERC).
Features	Siliconised Needle for less Trauma.
Type	2 part i.p. needle, 3 part i.p. needle
Size	18g*15 cms , 18g*20 cms
Use	For single use
Packing	Tyvek pouch as per EN-868-5- series standard
Sterilizations	With Ethylene Oxide EN-550 series standard

NMS-GW



GUIDE WIRE


Application	Guide Wires are used to gain ureteral access to establish a tract and to assist in the placement, replacement and exchange of medical devices during urological procedures
Features	Guide Wires are available in stainless steel PTFE coated , Nitinol hydrophilic coated, zebra and black & white.
Size	PTFE Guide wire : 0.018" / 0.025" / 0.032" / 0.035" / 0.038" Striped Zebra Guide Wire : 0.032" / 0.035" Nitinol Guide Wire : 0.035"
Type	Straight Tip / J-tip (PTFE Guide wire) Straight Tip (Striped Zebra Guide Wire / Nitinol Guide Wire)
Length	150cm
Use	For single use
Packing	Tyvek pouch as per EN-868-5- series standard
Sterilizations	With Ethylene Oxide EN-550 series standard

NMS-BDS



BALLOON DILATOR SET

Application	Used for ureteral dilation prior to ureteral stone manipulation or ureteroscopy, as well as for dilation of the intramural ureter.
Features	Shaft Size: 5 Fr - Balloon Dia: 5 mm - Balloon Length: 4 cm - Catheter Length: 75 cm - RBP: 14 Atm - Use Guide Wire: 0.035 inch - With Radiopaque markers
Use	For single use
Packing	Tyvek pouch as per EN-868-5- series standard
Sterilizations	With Ethylene Oxide EN-550 series standard



NIDHI[®] Meditech Systems
(Complete Urological Solutions Since-1995)

NIDHI CORPORATE HOUSE

No. H-28, H-29, H-30, Abhishek Estate-3,
Near Asarwa Over Bridge, Opp. Haripura Bus Stop, Haripura, Asarwa, Ahmedabad-380016
Cell / Whatsapp : +91-9510371980
Email : ndshah@nidhimeditech.com, info@nidhimeditech.com, sales@nidhimeditech.com
Website : www.nidhimeditech.com